




Ziehm Solo FD

Versatile design meets latest flat-panel technology

Image Quality: Changing the Game

The core of Medical Imaging: Image Quality

In medical imaging, the most important goal is obtaining the best clinical image quality. As the market leader in innovation, Ziehm Imaging provides a range of hardware and software features in our mobile C-arms, that helps to offer superior image quality. This changes the game for our clinicians and lets them discover new areas of clinical applications.



SmartDose

Best image quality. Minimized dose.

Comprehensive concept for dose reduction

Our latest improvements in SmartDose¹ help to display even the smallest details of complex anatomical areas and reduce dose with intelligent pulse regulation and optimized anatomical programs.

With significant dose savings, Ziehm Imaging sets the benchmark in user-friendly adjustment of dose exposure, and the SmartDose concept has been incorporated in the current generation of mobile C-arms.



Laser Positioning Device



Anatomical Programs



Low Dose Mode



Reduction of Pulse Frequency



High-Speed ADR



Removable Grid



Virtual Collimators



Automatic motion & position detection



Automatic adjustment for large patients



Exposure-free magnification



Beam Filtration?



ZAIP Algorithm & Filters

With its all-in-one design, the Ziehm Solo FD is one of the most compact C-arms on the market for even the smallest treatment rooms. The Ziehm Solo FD is equipped with the latest flat-panel technology, in either a 21 cm x 21 cm or a 31 cm x 31 cm detector. Perform a wide range of applications like orthopedics, trauma, pain management, peripheral vascular and anterior hips with excellent image quality. Versatile viewing options offer maximum flexibility and the latest improvements in SmartDose, the Ziehm Solo FD ensures best image quality at minimized dose.



Unique Selling Points

- Two detector sizes: 21 cm x 21 cm flat-panel or a 31 cm x 31 cm flat-panel
- 2.4 kW pulsed monoblock generator provides short, sharp pulses, producing razor-sharp images even if the patient is moving
- Simple positioning around patient and table due to the large C-arm opening
- Smallest footprint on the market even in limited room sizes
- Full size monitor mounted on C-arm
- Intuitive touchscreen user interface
- Comprehensive dose concept for high image quality and minimized dose

Imaging for a wide range of clinical applications

- Pain Management
- Orthopedics
- Peripheral Vascular
- Anterior Hips



IMAGING SOLUTIONS



Since 1972, Ziehm Imaging has been a company of firsts. Our visionary thinking has made us the mobile C-arm expert worldwide. We are committed to our mission of setting new technology standards. As an innovation leader, we were the first company to:

- Show Digital Subtraction Angiography (DSA) on a mobile C-arm
- Introduce a touchscreen user interface
- Incorporate a flat-panel detector on a mobile C-arm worldwide
- Launch a 3D C-arm with flat-panel detector on the market
- Introduce a mobile full-size C-arm with a CMOS detector

PROFESSIONAL TRAINING



Through our training programs, we empower our partners with the knowledge, skills, and insight they need to get the most out of our mobile C-arms. Our courses balance a classroom setting with hands-on experience to reveal the full benefits of our exceptional products.

Training is available at our Americas Academy Training facilities or customer's location (case-by-case basis).

NATIONWIDE SERVICE



The Customer Advantage:

- Nationwide coverage
- Fast, flexible & reliable service
- 24 hour delivery for parts inventoried in the US
- Comprehensive service covers the lifespan of your product

¹ In clinical practice, the use of SmartDose may reduce patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

² The technology Beam Filtration reduces dose exposure for Ziehm Imaging flat-detector systems in comparison with conventional filtration techniques. Data on File. Results may vary.

Ziehm Imaging
6280 Hazeltine National Dr | Orlando, FL 32822 | USA
Phone 800 503 4952 | Fax 321 445 5514
mail@ziehm.com | www.ziehm.com



Nationwide toll-free
**SERVICE & TECHNICAL
PHONE SUPPORT**
866.949.4346

